

# **SMAP: A Pipeline for Sample Matching in Proteogenomics**

**Version 1.0.0**

March 2021

# Contents

1. Introduction .....	3
1.1 Software requirement.....	3
1.2 Contact information .....	3
1.3 License .....	3
2. How to run SMAP (standard alone version) .....	4
2.1 Download the pipeline.....	4
2.2 Run SMAP program .....	4
2.3 Input data .....	4
2.3.1 Variant peptide table.....	4
2.3.2 Genotype file .....	5
2.3.3 Output files .....	5
3 Cloud-based SMAP .....	7
3.1 Introduction .....	7
3.2 Input data .....	7
3.3 Data summary & filtering .....	8
3.4 Run SMAP .....	10
4 References.....	11

## **1. Introduction**

SMAP is a pipeline designed for verifying and correcting sample identity for a large mass spectrometry (MS)-based proteomics project. SMAP takes a variant peptide data that can be generated using the proteogenomics approach. The program then infers allelic information for each sample based on its expression level of the variant peptides. The program finally aligns the MS-based proteomic samples with genomic information (i.e., genotypic data) by using two discriminant scores.

### **1.1 Software requirement**

SMAP has both standard alone and cloud-based versions. The standard alone version supports all 64-bit operating systems. The program is written by a combination of Perl and R. The minimum required Perl version should be Perl 5.6 or R 3.1.0.

### **1.2 Contact information**

For any questions, please contact Xusheng Wang ([xusheng.wang@und.edu](mailto:xusheng.wang@und.edu))

### **1.3 License**

This program is free software. You can redistribute and/or modify it under the terms of the GNU General Public License as published by the Free Software Foundation.

This program is distributed in the hope that it will be useful, but WITHOUT ANY WARRANTY; without even the implied warranty of MERCHANTABILITY or FITNESS FOR A PARTICULAR PURPOSE. See the GNU General Public License for more details.

## 2. How to run SMAP (standard alone version)

### 2.1 Download the pipeline

The pipeline could be downloaded from <https://github.com/UND-Wanglab/SMAP>

### 2.2 Run SMAP program

After installing SMAP program, you can run the program using the following command.

```
perl SMAP.pl -vf variant_peptide_table[file] -g genotype[file] -o result[file]
```

<code>--variant_peptide, -vf</code>	(A file containing quantitative values of variant peptides; required)
<code>--genotype, -g</code>	(A genotype file used sample verification; required )
<code>--output, -o</code>	(An output filename; required)
<code>--plex, -p</code>	(Multiplex number of the isobaric labeling approach)
<code>--fold_change, -fc</code>	(Signal to Noise ratio (optional; default is 3))
<code>--noise_level, -nl</code>	(The upper threshold of a noise level)
<code>--version, -h</code>	(Print version)
<code>--help, -h</code>	(Print help)
<code>--licence, -l</code>	(Print licencejump -s (search))

### 2.3 Input data

#### 2.3.1 Variant peptide table

The variant peptide table uses the following format:

Column 1: Peptide ID

Column 2: Gene/Protein

Column 3: Peptide Spectrum Match (PSM)

Column 4: SNP ID \*\*MUST MATCH GENOTYPE SNP ID

Column 5-N: Sample Peptide Quantification (One column per sample)

### An example of the variant peptide table

Peptide	Gene	PSM	SNP	2015-1341	...	2016-965	Internal standard
VSNEEKVR	CAPZA1	b20_f39.15855.1.3	chr1:113162494:G:A	53788.04	...	83146.90	46477.36
HWQQFYFLSTR	FBXO2	b20_f36.35042.1.3	chr1:11710561:T:G	25447.82	...	15590.47	19626.55
SIEDLLR	PDE4DIP	b20_f22.28382.1.2	chr1:144877111:G:T	13161.86	...	10127.43	8410.05

### 2.3.2 Genotype file

SMAP also takes a genotype in VCF format.

#### An example of the genotype data

#CHROM	POS	ID	REF	ALT	QUAL	FILTER	INFO	FORMAT	2014-2194	2014-2195	2014-2196
1	949608	chr1:949608:G:A	G	A	.	.	PR	GT	0/1	0/1	0/1
1	2441358	chr1:2441358:T:C	T	C	.	.	PR	GT	0/0	0/0	0/0
10	115644040	chr10:115644040:G:A	G	A	.	.	.	GT	0/1	0/0	0/1

### 2.3.3 Output files

SMAP generates a final report and several intermediate results.

The final report contains four columns, including Sample ID, Inferred ID, CSore and DeltaCScore.

#### An example of the final report

Sample ID	Inferred ID	CSore	DeltaCScore
2015-1341	2015-1341	4.22	0.70
2015-737	2015-737	4.03	0.56
2015-804	2015-804	3.70	0.59
2015-42	2015-37	3.14	0.51
2015-1555	2015-1555	2.91	0.54
2015-244	2015-244	2.62	0.44
2015-735	2015-735	2.53	0.43
2014-2200	2015-857	2.52	0.48
2016-958	2016-958	1.39	0.03
2016-965	2016-965	1.27	0.03
Internal standard	2015-1339	1.71	0.00

In addition, the program also generates three intermediate files, including sample-specific genotypes and inferred genotypes.

### An example of sample-specific genotype

#CHROM	POS	ID	REF	ALT	QUAL	FILTER	INFO	FORMAT	2014-2194	2014-2195	2014-2196
1	949608	chr1:949608:G:A	G	A	.	.	PR	GT	H	H	H
1	2441358	chr1:2441358:T:C	T	C	.	.	PR	GT	T	T	T
10	115644040	chr10:115644040:G:A	G	A	.	.	.	GT	C	H	H

### An example of inferred genotypes


SNP	2015-1341	2015-737	2015-804	2015-42	2015-1555	2015-244	2015-735	2014-2200	2016-958	2016-965	Internal standard
chr11:75298468:A:C	A	C	A	A	A	A	A	A	A	A	A
chr5:140503474:C:G	H	H	H	H	H	C	H	H	C	C	C
chr19:40408821:C:G	C	C	C	C	C	H	G	C	C	C	C

### 3 Cloud-based SMAP

You can visit the site at: <https://smap.shinyapps.io/smap/>

#### 3.1 Introduction

Navigation through the webpage can be done by clicking on any of the four tabs at the left.



#### 3.2 Input data

User can upload data using “Browse” buttons in “**Variant Peptide Table**” and “**Genotype Data**” menus. The format of both files can be found in section 2.3.

- The cloud-based SMAP application accepts .vcf files with meta-information (including none).
- Variant peptide data will be converted into log2 scale.


The screenshot shows two side-by-side input data upload forms. The left form is titled "Upload Variant Peptide Table" and the right form is titled "Upload Genotype Data". Both forms have a "Data Type" section with radio buttons for "User Data" and "Example Data" (selected). The left form also has a "Separator" section with radio buttons for "Comma", "Semicolon", and "Tab" (selected), and an "Is Input Data Log2 Normalized?" section with radio buttons for "Yes" and "No" (selected). Both forms have a "Choose input file" section with a "Browse..." button and "No file selected" text. The right form also has a "Choose Genotype File" section with a "Browse..." button and "No file selected" text. Below the forms, there is a note: "Input file MUST follow this column specification:" followed by a table with columns: Peptide, Gene/Protein, PSM, SNP, Sample 1 ..., Sample N.

### 3.3 Data summary & filtering


For both variant peptide and genotype data, SMAP provides summary values and relevant distributions for the input files.

- Large genotype files may take a few moments to load.

#### Variant peptide data summary (Example)



#### Genotype data summary (Example)





Users have the options to set each variant peptide filtering parameter (minimal intensity, maximal intensity, signal/noise ratio) and the genotype filtering parameter (number of missing genotypes tolerated per SNP) based on the data distributions.

### Variant Peptide Data Filtering

---

**Minimal Intensity Distribution (intensity of reference peptides)**

☐ No Filtering

☒ Input Filter:

Default value (if left blank above) is **17.45** (Mean + 2SD)

---

**Maximal Intensity Distribution**

☐ No Filtering

☒ Input Filter:

Default value (if left blank above) is **12.25** (Mean - 2SD)

---

**S/N Ratio Distribution**

☐ No Filtering

☒ Input Filter:

Default value (if left blank above) is **3**

---

### Genotype Data Filtering

---

**Filter Out SNPs With Missing Genotype Values**

☒ No Filter

☐ Filter

---

**Filter Out SNPs That DO NOT Appear in Filtered Variant Peptide Table**

☒ No Filter/Input VCF is Already Filtered

☐ Filter

***\*\*This filter will be applied when SMAP is run, but can be set here to see/download filtered vcf***

---

Default parameters are set (and selected if no user input is given) as follows:

- Minimal and maximal intensity filters are set based on the means and standard deviations of the minimal and maximal peptide distributions.
- Signal/noise ratio is always set at the default of 3.
- Number of missing genotypes tolerated filter is OFF at default, but the user can set this filter if they have a large amount of missing data.

### 3.4 Run SMAP

After selecting the desired filters, SMAP is by clicking the “Run SMAP” button.

SMAP

Introduction

Input Data

Data Summary & Filtering


Run SMAP

Run SMAP

Before hitting run: Ensure all filters are set to the desired values in previous tab

Run SMAP

After running SMAP, a table output will be generated displaying the variant peptide sample IDs and their matched genotype IDs. The Cscore and Delta Cscore for each match is also reported and graphed.



Users can download the results table by entering a desired file name and clicking download at the top right:

Download Results

Enter name for download file:

SMAP\_output

Download

## 4 References

1. Junmin Peng, J.E.E., Carson C Thoreen, Larry J Licklider, Steven P Gygi.(2003). Evaluation of multidimensional chromatography coupled with tandem mass spectrometry (LC/LC-MS/MS) for large-scale protein analysis the yeast proteome.pdf&gt;. J Proteome Res,2003(2),43-50.
2. Li, Y., Wang, X., Cho, J.H., Shaw, T.I., Wu, Z., Bai, B., Wang, H.,Zhou, S., Beach, T.G., Wu, G.\*, et al.\* (2016). JUMPg: An Integrative Proteogenomics Pipeline Identifying Unannotated Proteins in Human Brain and Cancer Cells. J Proteome Res,15(7), 2309-2320.
3. UniProt, C. (2021). UniProt: the universal protein knowledgebase in 2021. Nucleic acids research,49(D1), D480-D489.
4. Wang, K., Li, M., and Hakonarson, H. (2010). ANNOVAR: functional annotation of genetic variants from high-throughput sequencing data.Nucleic acids research,38(16), e164.
5. Wang, X., Li, Y., Wu, Z., Wang, H., Tan, H., and Peng, J. (2014). JUMP:a tag-based database search tool for peptide identification with high sensitivity and accuracy. Mol Cell Proteomics ,13(12), 3663-3673.